

Quantum Pathogen Scanner

Specifications

This technology was developed by NASA and the Russian space agency in order to find a device that can keep astronauts in good condition and protected against foreign pathogens during long trips. The machine works with a spectrum of radio waves (similar to the Rife wave generator) with the difference that it can also detect reflected waves to establish exactly at which wavelength foreign particles are reflected from. The machine can also determine the content both quantitatively (at which wavelength) and qualitatively (how many or large) the infestation is.

All known and “friendly” pathogens are programmed to be ignored and any content in the rest of the spectrum assumed foreign and potentially harmful.

Selection criteria of your Quantum Pathogen Scanner:

The conditions relating to your selection criteria are similar to that of the body scanner. At Art of Health we believe the extensive research for the correct product coupled with our superior technical support and commercial, product and application training forms the foundation of the added value you will be able to generate from your equipment.